

**Recombinant Mouse Nogo-66 Receptor/Reticulon 4 Receptor/NgR/RTN4R Protein(C-6His)**

**Cat.NO.: TP04830**

3th Edition

**Synonyms:**Reticulon-4 Receptor; Nogo Receptor; NgR; Nogo-66 Receptor; RTN4R; NOGOR

**Description:**Nogo Receptor (NgR) is a glycosylphosphoinositol (GPI)-anchored protein that belongs to the Nogo receptor family. Human NgR is predominantly expressed in neurons and their axons in the central nervous systems. As a receptor for myelin-derived proteins Nogo, myelin-associated glycoprotein (MAG) and myelin oligodendrocyte glycoprotein (OMG), NgR mediates axonal growth inhibition and may play a role in regulating axonal regeneration and plasticity in the adult central nervous system. NgR may be proposed as a potential drug target for treatment of various neurological conditions. Additionally, NgR may play a role in regulating the function of gap junctions.

**Form:**PBS

**Molecular Weight:**46.6 kDa

**Sequences:**Cys27-Ser447

**Purity:**> 95% by HPLC

**Concentration:**

**Endotoxin Level:**<1.0 EU per 1 ug of protein (determined by LAL method)

**Storage:**Can be stored at +4°C short term (1-2 weeks). For long term storage, aliquot and store at -20°C or -70°C. Avoid repeated freezing and thawing cycles.