

Recombinant Mouse Activin Receptor IB/Activin RIB/ALK-4/ACVR1B Protein(C-Fc)

Cat.NO.: TP04843

3th Edition

Synonyms:Activin Receptor Type-1B; Activin Receptor Type IB; ACTR-IB; Activin Receptor-Like Kinase 4; ALK-4; Serine/Threonine-Protein Kinase Receptor R2; SKR2; ACVR1B; ACVRLK4; ALK4

Description:Activin Receptor Type-1B (ACVR1B) is a single-pass type I membrane protein that belongs to the protein kinase superfamily. ACVR1B contains one GS domain and one protein kinase domain and is expressed in many tissues, most strongly in kidney, pancreas, brain, lung, and liver. ACVR1B acts as a transducer of activin or activin like ligands signals. Activin binds to either ACVR2A or ACVR2B and then forms a complex with ACVR1B, ACVR2A or ACVR2B activating ACVR1B through phosphorylation of its regulatory GS domain. They go on to recruit the R-SMADs, SMAD2 and SMAD3. ACVR1B also transducers signals of nodal, GDF-1, and Vg1. Mutations in ACVR1B are associated with pituitary tumors.

Form:PBS

Molecular Weight:37.7 kDa

Sequences:Leu32-Glu126

Purity:> 95% by HPLC

Concentration:

Endotoxin Level:<1.0 EU per 1 ug of protein (determined by LAL method)

Storage:Can be stored at +4°C short term (1-2 weeks). For long term storage, aliquot and store at -20°C or -70°C. Avoid repeated freezing and thawing cycles.