

Recombinant Human Parathyroid Hormone 1 Receptor/PTH1R Protein(Gly49, C-6His)

Cat.NO.: TP05048

3th Edition

Synonyms:Parathyroid hormone 1 receptor;PTH1R

Description:Parathyroid hormone 1 receptor(PTH1R) is a multi-pass membrane protein. The protein is expressed in high levels in bone and kidney and regulates calcium ionhomeostasis through activation of adenylate cyclase and phospholipase C. In bone, it is expressed on the surface of osteoblasts. When the receptor is activated through PTH binding, osteoblasts express RANKL (Receptor Activator of Nuclear Factor kB Ligand), which binds to RANK (Receptor Activator of Nuclear Factor kB) on osteoclasts. This turns on osteoclasts to ultimately increase the resorption rate.

Form:PBS

Molecular Weight:20.2 kDa

Sequences:Tyr23-Met 189

Purity:> 95% by HPLC

Concentration:

Endotoxin Level:<1.0 EU per 1 ug of protein (determined by LAL method)

Storage:Can be stored at +4°C short term (1-2 weeks). For long term storage, aliquot and store at -20°C or -70°C. Avoid repeated freezing and thawing cycles.