

Recombinant Human Ephrin-A3/EFNA3 Protein(C-6His)

Cat.NO.: TP05429

3th Edition

Synonyms:Ephrin-A3; EFL-2; EHK1 Ligand; EHK1-L; EPH-Related Receptor Tyrosine Kinase Ligand 3; LERK-3; EFNA3; EFL2; EPLG3; LERK3

Description:Ephrins-A3 belongs the Ephrins ligand family which involved in a variety of biological processes, especially in the nervous system and in erythropoiesis. It is shown that Ephrin-A3 is expressed in brain, skeletal muscle, spleen, thymus, prostate, testis, ovary, small intestine, and peripheral blood leukocytes. Ephrin-A3 has a GPI anchor following the extracellular sequence and a signal sequence of 22 amino acids. Ephrin-A3 can bind EphA2, EphA3, EphA4, EphA5, EphA6, EphA7, EphA8 and EphB1. Futhermore, it is associated with tumor growth and metastasis.

Form:PBS

Molecular Weight:22.3 kDa

Sequences:Gln23-Ser211

Purity:> 95% by HPLC

Concentration:

Endotoxin Level:<1.0 EU per 1 ug of protein (determined by LAL method)

Storage:Can be stored at +4°C short term (1-2 weeks). For long term storage, aliquot and store at -20°C or -70°C. Avoid repeated freezing and thawing cycles.