

Recombinant Human dUTP Pyrophosphatase/dUTPase Protein

Cat.NO.: TP05939

3th Edition

Synonyms:Deoxyuridine 5'-Triphosphate Nucleotidohydrolase Mitochondrial; dUTPase; dUTP Pyrophosphatase; DUT

Description:Deoxyuridine 5'-Triphosphate Nucleotidohydrolase Mitochondrial (dUTPase) belongs to the dUTPase family. dUTPase exists as a homotrimer and is involved in nucleotide metabolism. dUTPase produces dUMP, the immediate precursor of thymidine nucleotides and it decreases the intracellular concentration of dUTP so that uracil cannot be incorporated into DNA. The dUTPase increase in PCR product yield, length and fidelity enables further down-stream applications. These effects make dUTPase useful in PCR fidelity and yield-sensitive applications. dUTPase is specific for dUTP and is critical for the fidelity of DNA replication and repair.

Form:PBS

Molecular Weight:20.0 kDa

Sequences:Met 1-Asn164

Purity:> 95% by HPLC

Concentration:

Endotoxin Level:<1.0 EU per 1 ug of protein (determined by LAL method)

Storage:Can be stored at +4°C short term (1-2 weeks). For long term storage, aliquot and store at -20°C or -70°C. Avoid repeated freezing and thawing cycles.