



**Recombinant Human COL9A1 Protein (His Tag)**

**Cat.NO.: TP06377**

3th Edition

**Synonyms:**DJ149L1.1.2;EDM6;MED;STL4

**Description:**=VLOOKUP(A561,[1]Human!\$A:\$AG,33,FALSE)

**Form:**PBS

**Molecular Weight:**34.2 kDa

**Sequences:**Met 1-Pro328

**Purity:**> 95% by HPLC

**Concentration:**

**Endotoxin Level:**<1.0 EU per 1 ug of protein (determined by LAL method)

**Storage:**Can be stored at +4°C short term (1-2 weeks). For long term storage, aliquot and store at -20°C or -70°C. Avoid repeated freezing and thawing cycles.