

Recombinant Human FAM150A Protein (Fc Tag)

Cat.NO.: TP06809

3th Edition

Synonyms:FAM150A;UNQ9433

Description:

Form:PBS

Molecular Weight:15.0 kDa

Sequences:Met 1-Thr129

Purity:> 95% by HPLC

Concentration:

Endotoxin Level:<1.0 EU per 1 ug of protein (determined by LAL method)

Storage:Can be stored at +4°C short term (1-2 weeks). For long term storage, aliquot and store at -20°C or -70°C. Avoid repeated freezing and thawing cycles.