

**Recombinant Human DPEP2 Protein (Fc tag)**

**Cat.NO.: TP06848**

3th Edition

**Synonyms:**MBD2

**Description:**DPEP2 (MBD-2) belongs to the membrane-bound dipeptidase family. There are three members of this family as membrane-bound dipeptidase-1 (MBD-1),membrane-bound dipeptidase-2 (MBD-2) and membrane-bound dipeptidase-3 (MBD-3).MBD-2 is expressed at highest levels in lung, heart, and testis and at some what lower levels in spleen.MBD-2 is membrane-bound through a glycosylphosphatidyl-inositol linkage and probably is a metalloprotease which hydrolyzes leukotriene D4 (LTD4) into leukotriene E4 (LTE4).It is generally recognized that rapid cleavage of LTD4 is important in inactivating the broncho-and vaso-constrictive effects of cysteinyl LTs in asthmatic and inflammatory processes.

**Form:**PBS

**Molecular Weight:**65 kDa

**Sequences:**Met 1-Ser376

**Purity:**> 95% by HPLC

**Concentration:**

**Endotoxin Level:**<1.0 EU per 1 ug of protein (determined by LAL method)

**Storage:**Can be stored at +4°C short term (1-2 weeks). For long term storage, aliquot and store at -20°C or -70°C. Avoid repeated freezing and thawing cycles.