

Recombinant Human ERP72 / PDIA4 Protein (Fc tag)

Cat.NO.: TP06856

3th Edition

Synonyms:ERp-72;ERP70;ERP72

Description:ERP72, also known as PDIA4, is an endoplasmic reticulum luminal protein which belongs to the protein disulfide isomerase family. ERP72 is a stress protein and participates in the catalysis of protein-S-S-bond rearrangement. Both of PDIA4 and PDIA3 function as proteases, protein disulfide isomerases, phospholipases or an arrangement of these. ERP72 compose part of a large chaperone multiprotein complex comprising CABP1, DNAJB11, HSP90B1, HSPA5, HYOU, PDIA2, PDIA4, PPIB, SDF2L1, UGT1A1 and very small amounts of ERP29, but not, or at very low levels, CALR nor CANX.

Form:PBS

Molecular Weight:97.2 kDa

Sequences:Met 1-Thr641

Purity:> 95% by HPLC

Concentration:

Endotoxin Level:<1.0 EU per 1 ug of protein (determined by LAL method)

Storage:Can be stored at +4°C short term (1-2 weeks). For long term storage, aliquot and store at -20°C or -70°C. Avoid repeated freezing and thawing cycles.