

Recombinant Human FCGRT & B2M Heterodimer Protein

Cat.NO.: TP07099

3th Edition

Synonyms:alpha-chain;FCRN

Description:

Form:PBS

Molecular Weight:43.5 (31.8 + 11.7) kDa

Sequences:Met 1-Ser 297

Purity:> 95% by HPLC

Concentration:

Endotoxin Level:<1.0 EU per 1 ug of protein (determined by LAL method)

Storage:Can be stored at +4°C short term (1-2 weeks). For long term storage, aliquot and store at -20°C or -70°C. Avoid repeated freezing and thawing cycles.